








For scientific advice and support,
please contact us at :
info@u2bio.com



PROTEIN Section

Section	Recommended product	Cat No #	Size	Promotion
 Extraction	PRO-PREP™ Protein Extraction Solution (C/T)	17081	100 ml	4,100.-
	PRO-PREP™ Protein Extraction Solution (Cell/Tissue) (EDTA free)	17081.1	100 ml	4,100.-
 Measurement	PRO-MEASURE™ Protein Measurement Solution	21011	100 ml	4,200.-
Stacking Gel	0.5M Tris-Cl, pH 6.8 w/ SDS (Stacking gel buffer)	IBS-BT023	500 ml	1,350.-
Resolving Gel	1.5M Tris-Cl, pH 8.8 w/ SDS (Resolving gel buffer)	IBS-BT024	500 ml	1,530.-
Loading buffer	5X SDS-PAGE Loading Buffer w/2-Mercaptoethanol	IBS-BS024	10 ml	1,350.-
 Marker	GangNam Stain™ Prestained Protein Ladder (3 Colors)	24052	250 µl	2,800.-
Running	10X Tris-Glycine Buffer (w/ SDS)	IBS-BT014	1000 ml	2,160.-
Transfer	10x Western transfer buffer, carbonate	IBS-BW008	1000 ml	3,510.-
Blocking	Western blocking buffer	IBS-BW005	100 ml	1,530.-
 Stripping	NomeIRT™ Western Blot Stripping Buffer	21112	500 ml	3,500.-
Washing	10X TBS with Tween-20	IBS-BT007	1000 ml	3,060.-
 Detection	Miracle-Star™ Western Blot Detection Kit 🏆	16028	200 ml	4,200.-
	Fast-One Blocker (Blocking Solution and Signal Enhancer)	16174	500 ml	4,100.-



PRO-PREP™ Protein Extraction Solution (C/T) is;

Only takes 20min. for protein extraction from all kinds of cell and tissue and Easy to use: Single solution type kit

- Absorbance hindrance is not contained to eliminate error rate
- Protein stabilization buffer is included for protein stabilizing
- 5 kinds of protease inhibitor included to minimize protein degradation
- Monomeric formation is made by ionic detergent
- Useful for protein molecular weight measurement for western blot analysis



GangNam-STAIN™ Pre-stained Protein Ladder is;

- 12 bands between 10 and 245 kDa are available in a wide range of identifiable Pre-stained Protein Ladder
- Clear identification by two reference bands, 25 kDa(Green) and 75 kDa(Red)
- Easy-to-use products without prior processing
- Suitable for Western blotting, CBB staining, Silver staining and etc.
- High efficiency with small amount

PRO-MEASURE™ Protein Measurement Solution



- Bradford assay
- Simple application within 10 min

Miracle-Star™ Western Blot Detection Kit



NomeIRT™ Western Blot Stripping Buffer



- No smell and able to experiment at Clean hood
- Antibody can be removed by reaction at RT